

Buffers and Stocks

5X diluted Silica-coated Magnetic Beads

| Component | For 50 ml |
|------------------------------|-----------|
| Silica-coated Magnetic Beads | 10 ml |
| Add TE buffer to | 50 ml |

- [Home-Made Magnetic Beads](#)

TE buffer

| Component | Final conc. | For 50 ml |
|----------------------|-------------|-----------|
| 1 M Tris-HCl, pH 8.0 | 10 mM | 0.5 ml |
| 0.5 M EDTA, pH 8.0 | 1 mM | 0.1 ml |
| Add water to | | 50 ml |

1.5X GITC Lysis Buffer

| Component | Final conc. | For 50 ml |
|---------------------------------|-------------|---------------------|
| GITC (G54000) | 6 M | 35.46 g |
| 1 M Tris HCl, pH 8.0 | 75 mM | 3.75 ml |
| Sarkosyl | 3% | 1.5 g |
| EDTA | 30 mM | 3 ml of 0.5 M stock |

Adjust pH with HCl to 7.6-8.0 and adjust the volume with water to 50 ml

- Heat at 65 °C to dissolve.

TNES Buffer

| Component | Final conc. | For 50 ml |
|----------------------|-------------|-----------|
| 1 M Tris HCl, pH 8.0 | 100 mM | 5 ml |
| 5 M NaCl | 25 mM | 0.25 ml |
| 0.5 M EDTA | 10 mM | 1 ml |
| SDS | 10% w/v | 5 g |

- Add 25 µl of RNase A (10 mg/ml) per 1 ml TNES Buffer before use.

Wash Buffer

| Component | Final conc. | For 1 L |
|--------------|-------------|---------|
| Ethanol | 80% | 800 ml |
| Add water to | | 1 L |

Other materials

- Isopropanol

Procedures

1. Collect approximately 100 mg of plant tissue in two 2 ml collection tubes or 8 strip collection tube with two 4 mm beads. Freeze the samples at -80°C in a refrigerator or liquid nitrogen.

!!! info "For 8 strip collection tube, don't sample excessive tissue"

2. (**Highly recommend**) Lyophilize the samples overnight to enhance extraction, especially for root samples (Open the lid and cover with press 'n seal wrap and poke a hole in it).
3. Homogenize the samples using a MiniG at 1500 rpm for 45 seconds. Re-freeze the samples and repeat the homogenization process until all the samples are ground into a powder.

!!! info "Make sure the cap is completely sealed"

4. Add 200 μl of TNES Buffer with RNase A to the homogenized samples and mix by vortex.

!!! info "When using 8 strip collection tube, tilt-tap the tube first to prevent powder aggregate at the bottom"

5. Incubate the samples at 37°C for 15 minutes. A longer incubation time could degrade DNA.

!!! info "When using 8 strip collection tube, briefly spin down the sample to prevent cross contamination"

6. Add 400 μl of 1.5X GITC Lysis Buffer, mix by vortex, and incubate for 5 min at RT.

!!! info "When using 8 strip collection tube, briefly spin down the sample to prevent cross contamination"

7. Add 300 μl of Chloroform and mix by vortex.

8. Centrifuge at top speed for 5 min or 4000 g for 10 min (8 strip collection tube).

9. Prepare the following plates:

- **Sample-plate** (96-deep well plate{target="_blank"}):
 - 444 μL isopropanol
 - 222 μL 5X Silica-coated diluted magnetic beads (Mix well!!!!!!)
 - 333 μL DNA sample supernatant from step 8
- **IPA-plate** (96-deep well plate): 1 ml Isopropanol

- **Wash-plate-1** (96-deep well plate): 1 ml Wash Buffer
- **Wash-plate-2** (96-deep well plate): 1 ml Wash Buffer
- **Wash-plate-3** (96-deep well plate): 1 ml Wash Buffer
- **Elution-plate** ([Standard 96-well plate](#){target="_blank"}): 100 µL / 150 µl nuclease-free Water (after evaporation, only ~40 µl / 90 µl remain)
- **Tip-comb-plate:** [KingFisher 96 Tip Comb for DW Magnets](#){target="_blank"} in a Standard 96-well plate.

10. Run the **DNA_Extraction_General** protocol (1 hour) located in KingFisher Flex machine.

Notes

- In this method, magnetic beads were diluted with TE buffer, therefore, you can use multichannel pipette to dispense.