

Buffers and Solutions (Prepare freshly)

2 M NaOH

Component	Final conc.	For 100 mL
10 M NaOH	2 M	20 mL
Dissolve in		100 mL with water

0.1 M HCl

Component	Final conc.	For 60 mL
5 M HCl	0.1 M	1.2 mL
Dissolve in		60 mL with water

Fe₂/3 solution

Component	Final conc.	For 60 mL
FeCl₃ (VWR #470301-382) ^{target="_blank"}	0.333 M	3.24 g
FeCl₂ · 4 H₂O (VWR #470301-356) ^{target="_blank"}	0.167 M	2 g
Dissolve in	0.1 M HCl	60 mL

Other materials:

- Ethanol
- [Tetraethyl orthosilicate \(Sigma #86578-250ML\)](#)^{target="_blank"}
- [Ammonia solution \(Sigma #1054321011\)](#)^{target="_blank"}

Procedure

1. Degas 100 mL 2 M NaOH and 60 mL Fe₂/3 solution using a 0.22 µm filter.
2. Heat 100 mL 2 M NaOH in a 250 mL flat-bottom boiling flask with a stir bar at 100 rpm, 80°C using a magnetic stirrer.
3. Use a dripping funnel to add 60 mL Fe₂/3 solution dropwise into 100 mL 80°C 2 M NaOH with stirring at 400 rpm. Alternatively, slowly pour 60 mL Fe₂/3 solution into 100 mL 80°C 2 M NaOH along the edge.

!!! info "Note" * **Black precipitate of Fe₃O₄ (FeO·Fe₂O₃) is formed**

* ****Less oxygen, more yield****

4. Turn off heating and add 10 mL 25% ammonia solution, stir at 400 rpm for 30 min, and cool down to room temperature.
 5. Transfer the solution to a plastic container and settle the magnetic particles using a [strong magnet](#){target="_blank"}, discard the supernatant.
 6. Resuspend in 200 mL EtOH, pellet magnetically, and discard the supernatant. Repeat twice.
 7. Resuspend the pellet in 200 mL EtOH.
 8. Mix 1.8 L of 99% EtOH with 50 mL of 25% ammonia solution and 200 mL of the synthesized iron oxide MNPs (from step 7) in a heat-resistant 5 L bottle using a magnetic stirrer (400 rpm). Switch on the heating and allow the solution to heat up to 80 °C.
- !!! info "Note" * **Make sure the iron oxide MNPs evenly distributed in the solution, prevent large clusters**
9. Add 45 mL Tetraethyl orthosilicate (TEOS) under constant stirring and incubate for another 30 minutes.
 10. Add 400 mL of ddH₂O to the solution.
 11. Allow the reaction to proceed overnight.
 12. Cool down the solution to room temperature.
 13. Separate the coated MNPs using a strong neodymium magnet.
 14. Wash twice with pure water.
 15. Wash twice with pure ethanol.
 16. Resuspend the pellet with water to 50 mL.

!!! info "Note" * **Generate ~ 15-20 mL magnetic beads each reaction * The magnetic beads should be pure black**

Reference

1. <https://doi.org/10.1371/journal.pbio.3000107>
2. <https://bomb.bio/protocols/>